


## Analyze and Create Labeled video in Deeplabcut

1. Add the video to the deeplabcut project folder: videos
2. Open a conda environnement and type the following commands:  
**activate DEEPLABCUT**  
**ipython**  
**import deeplabcut**
3. To analyze the video, run the command: **deeplabcut.analyze\_video(...)**
4. To create a labeled video, run the command:  
**deeplabcut.create\_labeled\_video(...)**

## Convert the pixel data into cm using Espresso

For this step, each video orientation will be done separately.

1. Open Espresso
2. Upload the video in Espresso (File -> Video...)
3. Upload the config file (Model -> Config DLC)
4. Upload the corresponding csv file from the video folder in the deeplabcut project (Import -> CSV coordonnées XYs)
5. Click on the feature “calib” and place the scale on a known 10 cm object (round circle of the precision grip task)
6. Export the converted csv file (Export -> CSV coordonnées XY) while making sure to select the option of conversion from pixel to cm



## Calibration of cameras using Anipose

1. Take a calibration video with a charuco board. If needed, a script called [ChArUco\\_board\\_creation.py](#) can be used to create a charuco board and customize it's characteristics. If the number of squares or dictionary used to make the motif of the board is changed, the config.toml file must be changed.
2. Add the three calibration videos to the calibration folder of the anipose project and label them like such:

> This PC > OS (C:) > Users > Usagers > 2-Project1-Monika-2025-05-07 > ANIPOSE > 2024-11-05 > calibration

Name	Date	Type	Size	Length
calib-charuco-camA.avi	2025-08-07 2:12 PM	AVI Video File (...)	692,226 ...	00:01:54
calib-charuco-camB.avi	2025-08-06 10:18 AM	AVI Video File (...)	592,057 ...	00:01:54
calib-charuco-camC.avi	2025-08-07 2:12 PM	AVI Video File (...)	922,492 ...	00:01:54
calibration.toml	2025-08-11 4:00 PM	Toml Source File	2 KB	
detections.pickle	2025-08-11 3:59 PM	PICKLE File	8,733 KB	

3. Make sure the calibration.toml and detections.pickle file are deleted before running the calibration (they are the desired output of a calibration)
4. Change the date of the folder in which the calibration folder resides (in this example, 2024-11-05)
5. Activate a conda environment where anipose is installed and change the directory (cd) into the one with your project

```
(base) C:\Users\Usagers>activate test2
(test2) C:\Users\Usagers>cd 2-Project1-Monika-2025-05-07\ANIPOSE
```

6. Run the command: **anipose calibrate** in the conda prompt

```
(test2) C:\Users\Usagers\2-Project1-Monika-2025-05-07\ANIPOSE>anipose calibrate
Calibrating...
C:\Users\Usagers\2-Project1-Monika-2025-05-07\ANIPOSE\2024-11-05
C:\Users\Usagers\2-Project1-Monika-2025-05-07\ANIPOSE\2024-11-05\calibration\calibration.toml
C:\Users\Usagers\2-Project1-Monika-2025-05-07\ANIPOSE\2024-11-05\calibration\calib-charuco-camA.avi
100%|████████████████████████████████████████| 6894/6894 [03:56<00:00, 29.10it/s]
3070 boards detected
C:\Users\Usagers\2-Project1-Monika-2025-05-07\ANIPOSE\2024-11-05\calibration\calib-charuco-camB.avi
100%|████████████████████████████████████████| 6894/6894 [02:51<00:00, 40.30it/s]
2522 boards detected
C:\Users\Usagers\2-Project1-Monika-2025-05-07\ANIPOSE\2024-11-05\calibration\calib-charuco-camC.avi
100%|████████████████████████████████████████| 6894/6894 [05:57<00:00, 19.29it/s]
5925 boards detected
```

## Conversion of XY coordinates into XYZ using Anipose

1. In the folder pose-2d of the ANIPOSE folder, place the .h5 and the .pickle files generated by Deeplabcut for each video, as well as the .csv file for each video generated by Espresso. Change the names of the files so they correspond to the ones in the picture:

> This PC > OS (C:) > Users > Usagers > 2-Project1-Monika-2025-05-07 > ANIPOSE > 2024-11-05 > pose-2d

Name	Date modified	Type	Size
2024-11-05-vid01-camA.csv	2025-08-12 10:19 ...	Microsoft Excel ...	62,187 KB
2024-11-05-vid01-camA.h5	2025-07-14 1:42 P...	H5 File	65,246 KB
2024-11-05-vid01-camA.pickle	2025-07-14 1:42 P...	PICKLE File	2 KB
2024-11-05-vid01-camB.csv	2025-08-12 10:21 ...	Microsoft Excel ...	62,337 KB
2024-11-05-vid01-camB.h5	2025-08-11 2:39 P...	H5 File	65,246 KB
2024-11-05-vid01-camB.pickle	2025-08-11 2:39 P...	PICKLE File	2 KB
2024-11-05-vid01-camC.csv	2025-08-12 10:23 ...	Microsoft Excel ...	64,543 KB
2024-11-05-vid01-camC.h5	2025-08-11 4:22 P...	H5 File	65,246 KB
2024-11-05-vid01-camC.pickle	2025-08-11 4:22 P...	PICKLE File	2 KB

C1=camA, C2=camB, C3=camC

2. In the folder videos-raw in ANIPOSE, place the original, unlabeled videos. Change the names to correspond to the picture above.
3. In the conda prompt, run the command: **anipose triangulate**

```
(test2) C:\Users\Usagers\install-Monika\example_ANIPOSE>anipose triangulate
Triangulating points...
C:\Users\Usagers\install-Monika\example_ANIPOSE\2019-08-02\pose-3d\2019-08-02-vid01.csv
100%|██████████| 87900/87900 [01:20<00:00, 1092.83it/s]
C:\Users\Usagers\anaconda3\envs\test2\lib\site-packages\aniposelib\cameras.py:1253: NumbaWarning:
Compilation is falling back to object mode WITHOUT looplifting enabled because Function "_error_fun_triangular" failed type inference due to: No implementation
n of function Function(<function diff at 0x000002018B1DBBE0>) found for signature:

>>> diff(array(float64, 3d, C), n=int64, axis=Literal[int](0))

There are 2 candidate implementations:
- Of which 2 did not match due to:
  Overload in function 'np_diff_impl': File: numba\np\old_arraymath.py: Line 3538.
  With argument(s): '(array(float64, 3d, C), n=int64, axis=int64)':
  Rejected as the implementation raised a specific error:
    TypingError: got an unexpected keyword argument 'axis'
  raised from C:\Users\Usagers\anaconda3\envs\test2\lib\site-packages\numba\core\typing\templates.py:791

During: resolving callee type: Function(<function diff at 0x000002018B1DBBE0>)
During: typing of call at C:\Users\Usagers\anaconda3\envs\test2\lib\site-packages\aniposelib\cameras.py (1303)

File "C:\Users\Usagers\anaconda3\envs\test2\lib\site-packages\aniposelib\cameras.py", line 1303:
def _error_fun_triangular(self, params, p2ds,
    <source elided>
    # temporal constraint
    errors_smooth = np.diff(p3ds, n=n_deriv_smooth, axis=0).ravel() * scale_smooth
    ^
```

@jit(forceobj=True, parallel=True)					
Iteration	Total nfev	Cost	Cost reduction	Step norm	Optimality
0	1	2.6522e+15			9.97e+13
1	2	6.6763e+14	1.98e+15	1.28e+02	1.25e+13
2	3	1.6927e+14	4.98e+14	3.62e+04	1.55e+12
3	4	4.1686e+13	1.28e+14	2.08e+05	1.84e+11
4	5	8.1537e+12	3.35e+13	5.54e+05	1.93e+10
5	6	8.5552e+11	7.30e+12	7.20e+05	1.45e+09
6	7	6.3331e+10	7.92e+11	4.96e+05	8.74e+07

4. The output csv file with the XYZ coordinates will be in the newly created folder pose-3d.

## Parameters Analysis using VS Code

1. Set up a folder with the following files (the exact structure doesn't matter):
  - csv with XYZ,
  - csv with Espresso-converted data from the top camera,
  - a folder with the yaml files for that trial
2. Open Anaconda prompt of the dlcfix environment


```
(base) C:\Users\Usagers>activate dlcfix
(dlcfix) C:\Users\Usagers>cd C:\Users\Usagers\install-Monika\DL-Analysis\VS_Code\AnalyseDeeplabcut
```

3. Run the code you want by writing python (chosen\_code.py)

```
(dlcfix) C:\Users\Usagers\install-Monika\DL-Analysis\VS_Code\AnalyseDeeplabcut>python main.py
```

4. A user input window will appear. Press ok and select the right folders.

Instructions ×

 Select the folder of the project (1) and csv data (2) you would like to analyze.  
 Next input the labels that were used in the videos (4).  
 If you not, default paths and labels will be used.

OK

If running the main.py, use the csv with the xyz coordinates retrieved from anipose. If running the wrist\_distance.py, use the csv with the xy coordinates from the top camera. If running the success\_rate.py, the csv files used do not matter.

The description of each code's can be found in the Readme in the github.